

List of Validated Methods

Sl. No	Method Number	Molecule	Extraction Method	Calibration Range	Effective Date	Internal Standard	Anti coagulant	Instrument
1	STP006-02	Quetiapine	SPE	5 ng/mL to 1000 ng/mL	30 Oct 2010	Clozapine	K2 EDTA	LC-MS/MS 3200
2	STP008-01	Rabeprazole	LLE	5.068 ng/mL to 1000.651 ng/mL	01 Aug 2011	Quetiapine Hemifumarate	K2 EDTA	LC-MS/MS 3200
3	STP012-00	Metformin	SPE	25.0 ng/mL to 3500.0 ng/mL	04 Jun 2010	Phenformin	K2 EDTA	LC-MS/MS 3200
4	STP015-02	Cetirizine	SPE	0.500ng/mL to 100.124ng/mL	05 Jan 2011	Donepezil HCl	K2 EDTA	LC-MS/MS 3200
5	STP016-02	Omeprazole	LLE	4.999 ng/mL to 801.435 ng/mL	07 Jul 2011	Propranolol Hydrochloride	K2 EDTA	LC-MS/MS 3200
6	STP019-02	Nifedipine	LLE	0.500 ng/mL to 150.057 ng/mL	16 Jul 2011	Quetiapine Hemifumarate	K2 EDTA	LC-MS/MS 4000
7	STP021-01	Aprepitant	LLE	10.004 ng/mL to 5001.952 ng/mL	30 Jul 2011	Quetiapine Hemifumarate	K2 EDTA	LC-MS/MS 4000
8	STP023-01	Levetiracetam	PPE	0.255 µg/mL to 50.121 µg/mL	07 Feb 2012	Rabeprazole	K2 EDTA	LC-MS/MS 3200
9	STP026-01	Hydroxyzine and Cetirizine	SPE	0.250 ng/mL to 100.093 ng/mL & 1.000ng/mL to 501.179 ng/ml	28 Sep 2011	Quetiapine Hemifumarate	K2 EDTA	LC-MS/MS 4000
10	STP028-00	Bromocriptine	SPE	10.015 pg/mL to 513.322 pg/mL	28 Dec 2011	Quetiapine Hemifumarate	K2 EDTA	LC-MS/MS 3200
11	STP029-01	Montelukast	LLE	5.000ng/mL to 501.161ng/mL	15 Oct 2011	Rabeprazole	K2 EDTA	LC-MS/MS 3200
12	STP030-01	Phenelzine	SPE	505.075 pg/mL to 25213.521 pg/mL	19 Feb 2013	Hydroxyzine dihydrochloride	K2 EDTA	LC-MS/MS 3200
13	STP031-00	Lidocaine	LLE	0.303 ng/mL to 300.374 ng/mL	26 Jun 2012	Quetiapine Hemifumarate	K2 EDTA	LC-MS/MS 3200
14	STP034-00	Desvenlafaxine	LLE	401.580 ng/mL to 1.000 ng/mL	07 Apr 2012	Desvenlafaxine-D6	K2 EDTA	LC-MS/MS 3200
15	STP032-01	Levetiracetam	LLE	0.253 µg/mL to 35.060 µg/mL	24 Feb 2012	Trazadone	K2 EDTA	LC-MS/MS 4000
16	STP033-02	Scopolamine	SPE	20.800 pg/mL to 503.893 pg/mL	22 May 2013	Scopolamine D3	K2 EDTA	LC-MS/MS 3200
17	STP034-00	Desvenlafaxine	LLE	1.000 ng/mL to 401.580 ng/mL	10 Apr 2012	Desvenlafaxine-D6	K2 EDTA	LC-MS/MS 4000
18	STP035-01	Desvenlafaxine	LLE	1.002 ng/mL to 1000.165 ng/mL	21 May 2012	Desvenlafaxine-D6	K2 EDTA	LC-MS/MS 3200
19	STP037-01	Nitrofurantoin	LLE	10.010 ng/mL to 1502.225 ng/mL	03 Jun 2012	Fenofibric acid	K2 EDTA	LC-MS/MS 3200
20	STP039-00	Bupropion and Hydroxy Bupropion	SPE	1.004 ng/mL to 302.744 ng/mL and 1.017 ng/mL to 604.323 ng/mL	08 Aug 2012	Bupropion-D9	K2 EDTA	LC-MS/MS 3200
21	STP040-00	Diclofenac	SPE	100.001pg/mL to 12024.179pg/mL	23 Jun 2014	Diclofenac D4	K2 EDTA	LC-MS/MS 4000
22	STP041-02	Naproxen	LLE	0.502 µg/mL to 100.036 µg/mL	29 Dec 2012	Naproxen-D3	K2 EDTA	LC-MS/MS 4000
23	STP042-01	Ziprasidone	LLE	1.000 ng/mL to 350.026ng/mL	30 Sep 2012	Risperidone	K2 EDTA	LC-MS/MS 4000
24	STP044-02	Tetrabenazine and 2, 3-Dihydro Tetrabenazine	LLE	30.071pg/mL to 1202.782pg/mL and 2.008ng/mL to 120.177ng/mL	22 May 2013	Tetrabenazine-D7	K2 EDTA	LC-MS/MS 3200
25	STP045-00	Didanosine	SPE	10.0222g/mL to 3003.561 ng/mL	04 Sep 2012	Zidovudine	K2 EDTA	LC-MS/MS 4000
26	STP046-01	Fenofibric Acid	LLE	0.010µg/mL to 30.194 µg/mL and 0.010 µg/mL to 10.074 µg/mL	28 Nov 2012	Atorvastatin Calcium	K2 EDTA	LC-MS/MS 3200
27	STP048-01	Trientine and N1-Acetyl Trientine	SPE	10.007ng/mL to 3001.614ng/mL & 10.020 ng/mL to 3000.331ng/mL	28 Jul 2015	Trientine-D4 , N1-Acetyl Trientine-D4	K2 EDTA	LC-MS/MS 4000
28	STP049-00	Efavirenz, Emtricitabine & Tenofovir	SPE	50ng/mL to 5000ng/mL, 50ng/mL to 3500ng/mL and 10ng/mL to 800ng/mL.	17 Oct 2012	Zidovudine	K2 EDTA	LC-MS/MS 3200
29	STP051-01	Mesalamine	SPE	5.002 ng/mL to 700.375 ng/mL	04 Sep 2014	Mesalamine-D3	K2 EDTA	LC-MS/MS 4000
30	STP053-01	Methylphenidate	LLE	100.820pg/mL to 30048.200pg/mL	03 Jul 2013	Methylphenidate D9	K2 EDTA	LC-MS/MS 4000
31	STP054-01	Sulthiame	LLE	50.280ng/mL to 15001.920ng/mL	04 May 2013	Acetazolamide	K2 EDTA	LC-MS/MS 3200
32	STP055-00	Valsartan	LLE	20.008ng/mL to 10005.012 ng/mL	25 Mar 2013	Candesartan	K2 EDTA	LC-MS/MS 4000
33	STP056-00	Bupropion	LLE	1.010ng/mL to 501.760ng/mL and 2.010ng/mL to 1001.760ng/mL and 1.010 ng/mL to 500.440ng/mL and 1.010ng/mL to 500.520ng/mL.	28 May 2013	Bupropion d9	K2 EDTA	LC-MS/MS 4000
34	STP057-01	Raloxifene, Raloxifene 4-Glucuronide, Raloxifene 6-Glucuronide	SPE	20.005pg/mL to 1201.612pg/mL, 1.005ng/mL to 600.540ng/mL and 0.503ng/mL to 300.020ng/mL.	23 Nov 2013	Raloxifene D4,Raloxifene D4-6-Glucuronide	K2 EDTA	LC-MS/MS 4000
35	STP058-02	Etodolac	LLE	50.010ng/mL to 49992.520ng/mL	13 Mar 2014	Etodolac-D4	K2 EDTA	LC-MS/MS 4000
36	STP059-01	Carbidopa	SPE	2.004ng/mL to 400.847ng/mL	13 Sep 2013	Carbidopa D3	K2 EDTA	LC-MS/MS 3200
37	STP060-01	Clonidine	LLE	25.010pg/mL to 2001.600pg/mL	10 Oct 2013	Clonidine D4	K2 EDTA	LC-MS/MS 4000
38	STP061-00	Carbidopa and Levodopa	SPE	2.000 ng/mL to 400.010 ng/mL and 5.008 ng/mL to 2500.203 ng/mL	04 Nov 2013	Carbidopa D3,LevodopaD3	K2 EDTA	LC-MS/MS 3200

Sl. No	Method Number	Molecule	Extraction Method	Calibration Range	Effective Date	Internal Standard	Anti coagulant	Instrument
39	STP062-00	Flecainide	LLE	1.010ng/mL to 500.086ng/mL	03 Dec 2013	Flecainide Acetate D4	K2 EDTA	LC-MS/MS 3200
40	STP063-00	Lopinavir and Ritonavir	LLE	5.004 ng/mL to 2000.634 ng/mL and 1.005 ng/mL to 200.281 ng/mL	06 Nov 2013	Abacavir	K2 EDTA	LC-MS/MS 3200
41	STP064-00	Cefuroxime	LLE	50.013ng/mL to 10006.620ng/mL	03 Dec 2013	Cefotaxime	K2 EDTA	LC-MS/MS 4000
42	STP065-00	Nicotinic Acid and in Nicotinuric Acid	LLE	40.028ng/mL to 20006.28ng/mL and 20.031ng/mL to 10025.763ng/mL	27 Feb 2014	Nicotinic Acid-d4, Nicotinuric Acid-d4	K2 EDTA	LC-MS/MS 4000
43	STP066-01	Metaxalone	LLE	25.006ng/mL to 12059.526ng/mL	02 Mar 2014	Metaxalone-d6	K2 EDTA	LC-MS/MS 4000
44	STP067-01	Amlodipine	LLE	49.972pg/mL to 11999.280pg/mL	11 Apr 2014	Amlodipine-d4	K2 EDTA	LC-MS/MS 4000
45	STP068-01	Tetrabenazine and 2,3-Dihydro Tetrabenazine	LLE	10.007pg/mL to 1599.372pg/mL and 0.201ng/mL to 100.138ng/mL	03 May 2014	Tetrabenazine-D7,Alpha Hydroxy Tetrabenazine D7	K2 EDTA	LC-MS/MS 4000
46	STP069-01	Omeprazole	LLE	10.152ng/mL to 5000.098ng/mL	23 May 2014	Omeprazole-d3	K2 EDTA	LC-MS/MS 3200
47	STP070-01	Nevirapine	LLE	10.016ng/mL to 6002.551ng/mL	03 Jun 2014	Nevirapine-d3	K2 EDTA	LC-MS/MS 4000
48	STP072-01	Rivastigmine	LLE	0.201ng/mL to 80.354ng/mL	24 Jan 2015	Rivastigmine D6	K2 EDTA	LC-MS/MS 4000
49	STP076-00	Trientine and N1-Acetyl Trientine	SPE	10.006ng/mL to 3003.344ng/mL&10.017 ng/mL to 3000.777ng/mL	06 Jun 2015	Trientine-D4 , N1-Acetyl Trientine-D4	K2 EDTA	LC-MS/MS 3200
50	STP077-01	Nicotinic Acid and Nicotinuric Acid	LLE	10.005ng/mL to 5000.746ng/mL and 10.005ng/mL to 5001.985ng/mL	11 Dec 2014	Nicotinic Acid-d4,Nicotinuric Acid-d4	K2 EDTA	LC-MS/MS 3200/4000
51	STP080-01	Valsartan	LLE	20.004ng/mL to 15019.848ng/mL	31 Jan 2015	Valsartan-d3	K2 EDTA	LC-MS/MS 4000
52	STP081-01	Dofetilide	LLE	50.011pg/mL to 5000.874pg/mL	31 Jul 2015	Dofetilide D4	K2 EDTA	LC-MS/MS 4000
53	STP082-01	Losartan and Losartan Carboxylic Acid	LLE	5.002ng/mL to 1503.663ng/mL and 5.001ng/mL to 1500.881ng/mL	27 May 2015	Losartan-d4, Losartan Carboxylic Acid-d4	K2 EDTA	LC-MS/MS 4000
54	STP083-01	Buprenorphine	LLE	20.001pg/mL to 1600.957pg/mL	15 Sep 2015	Buprenorphine-D4	K2 EDTA	LC-MS/MS 3200
55	STP084-01	Diclofenac	SPE	100.033pg/mL to 12019.034pg/mL	13 May 2015	Diclofenac D4	K2 EDTA	LC-MS/MS 4000
56	STP085-01	Mirabegron	SPE	0.202ng/mL to 150.127ng/mL	30 Jul 2015	Mirabegron D5	K2 EDTA	LC-MS/MS 5500
57	STP086-01	Zolpidem	LLE	0.510ng/mL to 300.475ng/mL	29 Dec 2015	Zolpidem-D6	K2 EDTA	LC-MS/MS 4000
58	STP087-01	Mono Methyl Fumarate	SPE	40.000 ng/mL to 5001.621 ng/mL	22 Mar 2016	Mono Methyl Fumarate-D3	K2 EDTA	LC-MS/MS 4000
59	STP089-01	Clobazam and N-Desmethyl Clobazam	LLE	2.000ng/mL to 1200.000ng/mL and 1.000ng/mL to 600.000ng/mL	25 Aug 2015	Clobazam 13C ₆ , N Hydroxy Clobazam 13C ₆	K2 EDTA	LC-MS/MS 4000
60	STP090-01	Atovaquone	LLE	100.031 ng/mL to 50058.446 ng/mL	08 Nov 2016	Buparvaquone	K2 EDTA	LC-MS/MS 4000
61	STP092-01	Oxybutynin and N-Desethyl Oxybutynin	LLE	0.202 ng / mL to 50.262 ng / mL and 0.501 ng / mL to 250.083 ng / mL	03 Oct 2015	Oxybutynin-d11 and N-Desethyl oxybutynin-d5	K2 EDTA	LC-MS/MS 4000
62	STP093-01	Ranitidine	LLE	5.000ng/mL to 2000.970ng/mL	08 Feb 2016	Ranitidine-D6	K2 EDTA	LC-MS/MS 3200
63	STP094-01	Budesonide	SPE	20.006 pg/mL to 10019.268 pg/mL	05 Feb 2016	Budesonide-D8	K2 EDTA	LC-MS/MS 4000
64	STP095-01	Tetrabenazine and 2,3-DihydroTetrabenazine	LLE	10.002pg/mL to 1500.787pg/mL and 0.201ng/mL to 100.157ng/mL	18 Mar 2016	Tetrabenazine-D6, 2,3 DihydroTetrabenazine -D7	K2 EDTA	LC-MS/MS 4000
65	STP096-02	Desmopressin	SPE	2.001pg/mL to 500.526pg/mL	03 Oct 2016	Desmopressin-D5	K2EDTA	LC-MS/MS 4000
66	STP097-00	Azathioprine and 6-Mercaptopurine	SPE	0.500 ng/mL to 100.888 ng/mL and 1.000 ng/mL to 200.997ng/mL	18 Apr 2016	Azathioprine-D3,6-Mercaptopurine-13C2	K2EDTA	LC-MS/MS 4000
67	STP098-01	V-50(2, 2'Azobis (2-methylpropionamide) dihydrochloride)	NA	5.010 ng/mL to 1001.989 ng/mL	16 Jan 2016	N1-Acetyl Trientine	K2EDTA	LC-MS/MS 3200
68	STP099-00	Penicillamine	PPT	10.009 ng/mL to 5003.424 ng/mL	29 Mar 2016	Penicillamine-D3	K2EDTA	LC-MS/MS 4000
69	STP100-01	Fentanyl	LLE	5.003pg/mL to 2001.996pg/mL	06 Apr 2016	Fentanyl-D5	K2EDTA	LC-MS/MS 3200
70	STP101-01	Olmesartan	LLE	10.000 ng/mL to 3000.775 ng/mL	03 Jul 2016	Olmesartan-D6	K2EDTA	LC-MS/MS 4000
71	STP102-01	UDCA,TUDCA,GUDCA	SPE	50.034 ng/mL to 10002.455 ng/mL, 50.009 ng/mL to 10003.581 ng/mL and 2.002 ng/mL to 400.383 ng/mL	26 Dec 2016	Ursodiol-D4, Glycoursodeoxycholic-D4 Acid, Tauroursodeoxycholic-D4	K2EDTA	LC-MS/MS 4000
72	STP103-01	Glipizide	LLE	2.006 ng/mL to 1200.026 ng/mL	01 Nov 2016	Glipizide-D11	K3 EDTA	LC-MS/MS 3200
73	STP104-01	Sulfamethoxazole and Trimethoprim	LLE	500.100ng/mL to 50038.377 ng/mL and 50.020ng/mL to 4998.743ng/mL	05 Oct 2016	Sulfamethoxazole-D4,Trimethoprim-D3	K3 EDTA	LC-MS/MS 3200
74	STP105-01	Molindone	LLE	0.103 ng/mL to 100. 461 ng/mL	21 Oct 2016	Molindone -D8	K3 EDTA	LC-MS/MS 3200
75	STP106-01	Darifenacin	LLE	20.009 pg/mL to 20004.924 pg/mL	22 Oct 2016	Darifenacin -D4	K3 EDTA	LC-MS/MS 3200
76	STP107-02	Omeprazole Impurities	NA	0.202 ppm to 5.058 ppm & 0.202 ppm to 5.057 ppm	25 Sep 2016	NA	K3 EDTA	LC-MS/MS 3200
77	STP108-01	Desmopressin	SPE	2.001pg/mL to 500.000pg/mL	08 Nov 2016	Desmopressin-D5	K3 EDTA	LC-MS/MS 3200

Sl. No	Method Number	Molecule	Extraction Method	Calibration Range	Effective Date	Internal Standard	Anti coagulant	Instrument
78	STP109-01	Dexmethylphenidate	LLE	100.005 pg/mL to 50013.706 pg/mL	24 Oct 2016	Methylphenidate –D9	K3 EDTA	LC-MS/MS 4000
79	STP110-01	Methyl Prednisolone	LLE	5.001 ng/mL to 800.275 ng/mL	06 Apr 2017	Methyl Prednisolone –D3	K3EDTA	LC-MS/MS 3200
80	STP111-01	Rosuvastatin	LLE	0.501 ng/mL to 100.519 ng/mL	12 May 2017	Rosuvastatin-D6	Whole blood	LC-MS/MS 4000
81	STP112-01	Esomeprazole	LLE	5.002 ng/mL to 4000.755 ng/mL	06 Oct 2017	Omeprazole-D3	Whole Blood	LC-MS/MS 4000
82	STP113-03	Ranolazine	LLE	10.002 ng/mL to 5000.495 ng/mL	30 Apr 2018	Ranolazine –D8	K2 EDTA	LC-MS/MS 3200
83	STP114-02	Ezetimibe and Ezetimibe Phenoxy Glucuronide	SPE	0.203ng/mL to 10.009 ng/mL and 1.001ng/mL to 500.215ng/mL	12 Dec 2017	Ezetimibe-D4,Ezetimibe Phenoxy Glucuronide-D4	K2 EDTA	LC-MS/MS 4000
84	STP115-01	Prednisolone and Trimeprazine	LLE	0.200 ng/mL to 200.074 ng/mL and 0.021 ng/mL to 20.663 ng/mL	10 Oct 2018	Prednisolone-D6, Trimeprazine-D6	K2 EDTA	LC-MS/MS 3200
85	STP116-00	Trans-Phytonadione and Cis-Phytonadione	LLE	0.502 ng/mL to 400.727 ng/mL and 0.501 ng/mL to 80.266 ng/mL	21 Nov 2017	Trans Phytonadione-D7,Cis Phytonadione –D7	K2 EDTA	LC-MS/MS 3200
86	STP117-01	Tacrolimus	LLE	0.208 ng/mL to 99.996 ng/mL	14 Dec 2017	Sirolimus	K2 EDTA	LC-MS/MS 4000
87	STP118-01	Simvastatin and Simvastatin Hydroxy Acid	LLE	0.201 ng/mL to 80.215 ng/mL and 0.101 ng/mL to 40.047 ng/mL	02 Dec 2017	Simvastatin-D6 ,Simvastatin Hydroxy acid-D6	K2 EDTA	LC-MS/MS 3200
88	STP119-02	Amoxicillin and Clavulanic acid	SPE	400.053 ng/mL to 60011.595 ng/mL and 100.018 ng/mL to 15007.755 ng/mL	26 Jul 2017	Amoxicillin-D4,Cefotaxime	K3EDTA	LC-MS/MS 3200
89	STP120-01	Cyclosporine	LLE	10.006 ng/mL to 2000.833ng/mL	01 Oct 2018	Cyclosporine A D4	K2 EDTA	LC-MS/MS 4000
90	STP121-02	Celecoxib	LLE	10.011ng/mL to 3003.883 ng/mL	09 Nov 2017	Celecoxib –D4	K2 EDTA	LC-MS/MS 4000
91	STP122-01	Amoxicillin and Clavulanic acid	SPE	400.004 ng/mL to 60028.120 ng/mL and 100.026 ng/mL to 15004.429 ng/mL	26 Jul 2017	Amoxicillin-D4,Cefotaxime	K2 EDTA	LC-MS/MS 4000
92	STP123-01	Fluoxetine	LLE	0.501 ng/mL to 100.062 ng/mL	28 Jul 2017	Fluoxetine-D5	K2 EDTA	LC-MS/MS 3200/4000
93	STP124-01	Morphine and Morphine 6β-D-glucuronide	SPE	0.501 ng/mL to 200.000 ng/mL and 0.501 ng/mL to 200.000 ng/mL	19 Nov 2018	Morphine –D6 ,Morphine 6β-D-glucuronide -D3	NA	LC-MS/MS 3200
94	STP125-01	Itraconazole and Hydroxy Itraconazole	LLE	1.009 ng/mL to 500.520 ng/mL and 1.010 ng/mL to 500.677 ng/mL	18 Nov 2017	Itraconazole-D9,Hydroxy Itraconazole-D8	K2 EDTA	LC-MS/MS 4000
95	STP126-01	Sulfamethoxazole and Trimethoprim	LLE	500.189ng/mL to 50016.472 ng/mL and 10.019ng/mL to 1001.863ng/mL	31 Aug 2017	Sulfamethoxazole-D4,Trimethoprim-D3	Sodium Heparin	LC-MS/MS 5500
96	STP127-01	Deracoxib	LLE	5.005 ng/mL to 2500.730 ng/mL	21 Mar 2018	Deracoxib-D3	Sodium Heparin	LC-MS/MS 5500
97	STP128-00	Chlorpromazine and 7-Hydroxy Chlorpromazine	LLE	0.100 ng/mL to 50.105 ng/mL and 0.101 ng/mL to 50.396 ng/mL	19 Nov 2018	Chlorpromazine -D6,7-Hydroxy Chlorpromazine -D6	Sodium Heparin	LC-MS/MS 4500
98	STP129-01	Propranolol and 4-Hydroxy Propranolol	LLE	2.000 ng/mL to 400.212 ng/mL and 1.000 ng/mL to 200.080 ng/mL	29 Aug 2017	Propranolol -D7,4-Hydroxy Propranolol -D7	K2 EDTA	LC-MS/MS 5500
99	STP130-01	Meloxicam	LLE	5.033 ng/mL to 3002.616 ng/mL	12 Nov 2018	Meloxicam -D3	K2 EDTA	LC-MS/MS 8060
100	STP131-02	Mycophenolic acid	LLE	100.002 ng/mL to 40005.318ng/mL	23 Jan 2018	Mycophenolic acid –D3	K2 EDTA	LC-MS/MS 3200
101	STP132-00	Bile Acid	NA	0.050 mM to 29.965 mM	14 Nov 2018	NA	NA	BA/HPLC/002
102	STP133-02	Mycophenolate Mofetil	LLE	30.005 pg/mL to 30031.740 pg/mL	03 Apr 2018	Mycophenolate Mofetil-D4	K2 EDTA	LC-MS/MS 5500
103	STP134-02	Rufinamide	LLE	50.033 ng/mL to 10004.941 ng/mL	20 Oct 2017	Rufinamide –D2	K2 EDTA	LC-MS/MS 4000
104	STP135-01	Marbofloxacin	LLE	100.124 ng/mL to 12002.920 ng/mL	18 Nov 2017	Marbofloxacin –D5	K2 EDTA	LC-MS/MS 5500
105	STP136-01	Carprofen	LLE	0.502 µg/mL to 50.026 µg/mL	15 Dec 2017	Ibuprofen	K2 EDTA	LC-MS/MS 4000
106	STP137-01	Solifenacin	LLE	0.200 ng/mL to 100.675 ng/mL	11 Jan 2018	Solifenacin-D5	K2 EDTA	LC-MS/MS 3200
107	STP138-01	Marbofloxacin	LLE	100.135 ng/mL to 12001.383 ng/mL	28 Dec 2017	Marbofloxacin-D5	K2 EDTA	LC-MS/MS 8060
108	STP139-01	Abiraterone	LLE	0.200 ng/mL to 200.158 ng/mL	21 Dec 2017	Abiraterone-D4	K2 EDTA	LC-MS/MS 3200
109	STP141-01	Bumetanide	LLE	0.507 ng/mL to 152.035 ng/mL	01 Feb 2018	Bumetanide-D5	K2EDTA	LC-MS/MS 3200
110	STP142-01	Ethacrynic Acid	LLE	2.000 ng/mL to 400.042 ng/mL	02 Jan 2018	Ethacrynic Acid-D5	K2EDTA	LC-MS/MS 4000
111	STP143-01	Pimobendan	LLE	0.100 ng/mL to 20.019 ng/mL	26 Sep 2018	Pimobendan-13CD3	K2EDTA	LC-MS/MS 4000
112	STP144-01	Edetic Acid	NA	0.250 ppm to 5.000 ppm	04 Apr 2018	NA	Topical solution	LC-MS/MS 4000
113	STP145-01	Amino Acids	NA	0.020 µg/mL to 15.296 µg/mL	14 Nov 2018	NA	Glatiramer acetate	LC-MS/MS 4500
114	STP146-01	Paliperidone	LLE	0.100 ng/mL to 50.056 ng/mL	27 Feb 2018	Paliperidone-D4	K2EDTA	LC-MS/MS 4000
115	STP147-01	Bromocriptine	LLE	10.005 pg/mL to 1500.814 pg/mL	10 Apr 2018	Bromocriptine-13CD3	K2EDTA	LC-MS/MS 5500
116	STP148-02	Oxcarbazepine and 10,11-Dihydro-10-Hydroxy Carbamazepine	LLE	10.166 ng/mL to 5080.460 ng/mL and 55.300 ng/mL to 25086.130 ng/mL	06 Jul 2018	Oxcarbazepine- D4&10,11-Dihydro-10-Hydroxy Carbamazepine-D4	K2EDTA	LC-MS/MS 3200
117	STP149-00	Desfuroylceftiofur acetamide	PPT	0.101 µg/mL to 50.144 µg/mL	30 Apr 2018	Desfuroylceftiofur acetamide-D3	Lithium Heparin	LC-MS/MS 4000

Sl. No	Method Number	Molecule	Extraction Method	Calibration Range	Effective Date	Internal Standard	Anti coagulant	Instrument
118	STP150-01	Enrofloxacin	LLE	1.001 ng/mL to 500.962 ng/mL	09 Apr 2018	Enrofloxacin-D5	K2EDTA	LC-MS/MS 8060
119	STP151-01	Enrofloxacin	LLE	1.000 ng/mL to 500.329 ng/mL	27 Mar 2018	Enrofloxacin-D5	K2EDTA	LC-MS/MS 8060
120	STP154-01	Eplerenone	LLE	5.002 ng/mL to 2500.556 ng/mL	08 Jan 2018	Eplerenone-D3	K2EDTA	LC-MS/MS 4500
121	STD155-01	Teriflunomide	LLE	20.001 ng/mL to 3000.511 ng/mL	30 May 2018	Teriflunomide-D4	K2EDTA	LC-MS/MS 4000
122	STP156-00	Metformin	LLE	25.047 ng/mL to 2515.756 ng/mL	29 May 2018	Metformin-D6	K2EDTA	LC-MS/MS 4000
123	STD158-01	Droxidopa	SPE	10.001 ng/mL to 3000.036 ng/mL	03 Jul 2018	Droxidopa-D3	K2EDTA	LC-MS/MS 4500
124	STD159-00	Erythromycin	LLE	20.008 ng/mL to 4003.467 ng/mL	27 Aug 2018	Erythromycin -13C-D3	K2EDTA	LC-MS/MS 3200
125	STD160-01	Cyclophosphamide	LLE	0.100 µg /mL to 30.013 µg/mL	10 Aug 2018	Cyclophosphamide-D8	K2EDTA	LC-MS/MS 4000
126	STD161-01	Eliglustat	LLE	0.201 ng/mL to 100.485 ng/mL	18 Jul 2018	Eliglustat-D4	K2EDTA	LC-MS/MS 4500
127	STD163-00	Carbamazepine	LLE	50.007 ng/mL to 10006.700 ng/mL	01 Sep 2018	Carbamazepine-D10	K2EDTA	LC-MS/MS 3200
128	STD165-00	Erythromycin Ethyl Succinate	SPE	20.006 ng/mL to 4001.736 ng/mL	13 Aug 2018	Erythromycin Ethyl Succinate-D6	K2EDTA	LC-MS/MS 3200
129	STD166-01	Aminocaproic Acid	PPT	1.002 µg/mL to 200.172 µg/mL	12 Oct 2018	Aminocaproic Acid -D6	K2EDTA	LC-MS/MS 3200
130	STD167-00	Azithromycin	LLE	5.000 ng/mL to 1500.370 ng/mL	31 Aug 2018	Erythromycin	K2EDTA	LC-MS/MS 4000
131	STP168-00	Citric acid	PPT	5.016 µg / mL to 2019.568 µg / mL	08 Nov 2018	Citric acid-D4	Urine	LC-MS/MS 3200
132	STP169-01	Bisoprolol	LLE	0.510 ng/mL to 151.092 ng/mL	19 Sep 2018	Bisoprolol-D5	K2EDTA	LC-MS/MS 4500
133	STP170-00	Potassium	LLE	51.107 µg/mL to 3005.233 µg/mL	05 Nov 2018	Barium	Urine	ICP-MS NexION 2000
134	STP171-00	Carbidopa and Levodopa	LLE	2.009 ng/mL to 400.951 ng/mL and 5.015 ng/mL to 2501.636 ng/mL	15 Sep 2018	Carbidopa-D3 & Levodopa-D3	K2EDTA	LC-MS/MS 5500
135	STP172-01	Butalbital and Codeine	LLE	25.015 ng/mL to 5008.000 ng/mL and 1.018 ng/mL to 200.000 ng/mL	15 Oct 2018	Butalbital -D5 and Codeine -D6	K2EDTA	LC-MS/MS 4500